

# SETHNESS PRODUCTS COMPANY

## P60

### LIQUID CARAMEL COLOR

#### CLASS III – E150c

### SPECIFICATION

<b>Description</b>	According to Title 21 CFR 73.85, the color additive Caramel is the dark brown liquid or solid material resulting from the carefully controlled heat treatment of food-grade carbohydrates. Certain food-grade acids, alkalis and salts may be employed to assist caramelization. This product is GRAS (21 CFR 182.1235).		
<b>Ingredient Statement</b>	Caramel Color, Caramel, Caramel Color (100%) or Caramel (100%)		
<b>Organoleptic Properties</b>	Characteristic bitter, burnt sugar taste.		
<b>GMO Status</b>	This product is GMO-free by test (PCR negative). This Caramel Color is produced with genetically engineered components.		
<b>Chemical/Physical Properties (as manufactured)</b>	<b>Property</b>	<b>Range</b>	
	Tinctorial Power, $K_{0.56}$ (560nm)	0.055-0.070	
	Baume' @ 60°F	35.5-36.5	
	Specific Gravity @ 60°F	1.324-1.336	
	Density (lb./gal. @ 60°F)	11.03-11.13	
	pH "as is"	4.5-5.5	
	Color Intensity (typical) (610nm)	0.028-0.037	
<b>Microbiological Properties (tested on audit basis)</b>	<b>Test</b>	<b>Unit</b>	<b>Limit</b>
	Aerobic Plate Count (APC)	cfu/g, max.	200
	Yeast	cfu/g, max.	10
	Mold	cfu/g, max.	10
	E. coli	MPN/g	<3
	Salmonella	cfu/25 g	Negative
<b>Packaging</b>	4x1 gal./case, 5 gal. containers, 55 gal. drum and bulk		
<b>Storage Conditions</b>	Cool and dry environment, preferably not to exceed 90°F. Product should not be allowed to freeze.		
<b>Shelf-life</b>	Two year minimum		
<b>Good Manufacturing Practice</b>	Rotate inventory and maintain stock no longer than one year.		

Please see corresponding Nutritional document for additional information.

P60 SPEC 010/LAB/LIQ	PAGE 1 OF 1	Authorized by JME
1347 Beaver Channel Parkway, Clinton, IA 52732	Original: 01/15/04	Revised: 04/25/19

NOTE: Once this document is copied, printed, emailed, downloaded, or modified in any way from what is stated on its' docsign sheet, it becomes UNCONTROLLED. For the current Controlled Document or assistance contact [lab@sethness.com](mailto:lab@sethness.com) or [www.sethness.com](http://www.sethness.com) or call 888-772-1880.